

Lateral Flow Mini- Materials Kit (MATL-005)

Item	Application	Supplier	Quantity/Dimensions
CN95	Membrane	Sartorius	25 Strips, 300mm x 25mm
CN140	Membrane	Sartorius	25 Strips, 300mm x 25mm
FF120HP	Membrane	Cytiva	25 Strips, 300mm x 25mm
FF170HP	Membrane	Cytiva	25 Strips, 300mm x 25mm
6614	Conjugate Pad	Ahlstrom	2 Sheets, 200mm x 300mm
8951	Conjugate Pad	Ahlstrom	2 Sheets, 200mm x 300mm
8964	Sample Pad	Ahlstrom	2 Sheets, 200mm x 300mm
222	Absorbent Pad	Ahlstrom	2 Sheets, 200mm x 300mm
Backing Cards	Adhesive Backing Support	DCN Dx	50 pieces, 60mm x 300mm
Backing Cards	Adhesive Backing Support	DCN Dx	50 pieces, 80mm x 300mm
Silica Gel Pouch	Desiccator	DCN Dx	200 pieces, 0.5g/pouch
Foil Pouch	Foil Pouch	DCN Dx	100 pieces, 3" x 4.5"

Please note: Membranes may be substituted with similar-performance materials.. These substitutions are reviewed and approved by DCN Dx scientists to ensure comparable early stage performance, but may vary slightly in chemistry or flow characteristics.

Thank you for your interest in our DCNovations Lateral Flow Materials Kits!

Place your order online at www.dcn dx.com/store or contact us with any questions at customercare@dcndx.com