

Lateral Flow Materials Kit (MATL-025)

Item	Application	Supplier	Quantity/Dimensions
CN95	Membrane	Sartorius	25 Strips, 300mm x 25mm
CN140	Membrane	Sartorius	25 Strips, 300mm x 25mm
FF120HP	Membrane	Cytiva	25 Strips, 300mm x 25mm
FF170HP	Membrane	Cytiva	25 Strips, 300mm x 25mm
FF80HP	Membrane	Cytiva	25 Strips, 300mm x 25mm
FF120HP Plus	Membrane	Cytiva	25 Strips, 300mm x 25mm
FF80HP Plus	Membrane	Cytiva	25 Strips, 300mm x 25mm
FF170HP Plus	Membrane	Cytiva	25 Strips, 300mm x 25mm
6614	Conjugate Pad	Ahlstrom	2 Sheets, 200mm x 300mm
8951	Conjugate Pad	Ahlstrom	2 Sheets, 200mm x 300mm
6615	Conjugate Pad	Ahlstrom	2 Sheets, 200mm x 300mm
1660	Sample Pad	Ahlstrom	2 Sheets, 200mm x 300mm
8964	Sample Pad	Ahlstrom	2 Sheets, 200mm x 300mm
Standard 17	Sample Pad	Cytiva	25 Strips, 300mm x 22mm
Fusion 5	Sample Pad	Cytiva	25 Strips, 300mm x 25mm
222	Absorbent Pad	Ahlstrom	2 Sheets, 200mm x 300mm
320	Absorbent Pad	Ahlstrom	2 Sheets, 200mm x 300mm
319	Absorbent Pad	Ahlstrom	2 Sheets, 200mm x 300mm
Backing Cards	Adhesive Backing Support	DCN Dx	100 pieces, 60mm x 300mm
Backing Cards	Adhesive Backing Support	DCN Dx	100 pieces, 80mm x 300mm
Plastic Cassettes	Test Strip Housing	DCN Dx	100 pieces, suitable for 5mm test strip
Silica Gel Pouch	Desiccator	DCN Dx	200 pieces, 0.5g/pouch
Foil Pouch	Foil Pouch	DCN Dx	100 pieces, 3" x 4.5"

Please note: Membranes may be substituted with similar-performance materials.. These substitutions are reviewed and approved by DCN Dx scientists to ensure comparable early stage performance, but may vary slightly in chemistry or flow characteristics.